

MICROPIPETTES



- Fixed or variable volume, general purpose micropipettes for sampling and dispensing precise amount of liquid in μl volumes
- Operates on the air displacement principle and uses disposable microtips
- Cover the volume range of $0.5 \mu\text{l}$ to $1000 \mu\text{l}$
- All models are equipped with a built in tip ejector



ORDERING INFORMATION FOR MICROPIPETTES

Cat. No.	Product	Volume
FPP5	BhatBiopette - 'Premium' Fixed Micropipette	$5 \mu\text{l}$
FPP10	BhatBiopette - 'Premium' Fixed Micropipette	$10 \mu\text{l}$
FPP20	BhatBiopette - 'Premium' Fixed Micropipette	$20 \mu\text{l}$
FPP50	BhatBiopette - 'Premium' Fixed Micropipette	$50 \mu\text{l}$
FPP100	BhatBiopette - 'Premium' Fixed Micropipette	$100 \mu\text{l}$
FPP500	BhatBiopette - 'Premium' Fixed Micropipette	$500 \mu\text{l}$
FPP1000	BhatBiopette - 'Premium' Fixed Micropipette	$1000 \mu\text{l}$
VPP10	BhatBiopette - 'Premium' Variable Micropipette	$0.5 - 10 \mu\text{l}$
VPP20	BhatBiopette - 'Premium' Variable Micropipette	$2 - 20 \mu\text{l}$
VPP50	BhatBiopette - 'Premium' Variable Micropipette	$5 - 50 \mu\text{l}$
VPP100	BhatBiopette - 'Premium' Variable Micropipette	$10 - 100 \mu\text{l}$
VPP200	BhatBiopette - 'Premium' Variable Micropipette	$20 - 200 \mu\text{l}$
VPP1000	BhatBiopette - 'Premium' Variable Micropipette	$100 - 1000 \mu\text{l}$